



PROVOX[®] HME Life

Figures

Disclaimer

Atos Medical offers no warranty – neither expressed nor implied – to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patent and trademarks

Provox[®] is a registered trademark owned by Atos Medical AB, Sweden.

Provox[®] Life[™] is a trademark of Atos Medical AB.

For information about protective rights (e.g. patents), please refer to our webpage www.atosmedical.com/patents

Manufacturer; Hersteller; Fabrikant; Fabricant; Produttore;
 Fabricante; Fabricante; Tillverkare; Producent; Produsent;
 Valmistaja; Framleiðandi; Tootja; Ražotājs; Gamintojas;
 Výrobce; Gyártó; Výrobca; Proizvajalec; Producent; Producător;
 Proizvođač; Proizvođač; Κατασκευαστής; Производител; Üretici;
 მწარმოებელი; Ἀρτοποιητής; İstehsalçı; Производитель; Pro-
 dusen; Pengilang; 製造業者; 제조사; 製造商; 制造商;

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1. Descriptive Information

1.1 Intended use

Provox Life™ HMEs are single-use heat and moisture exchangers for patients breathing through a tracheostoma.

1.2 CONTRAINDICATIONS

The product shall not be used by patients with a decreased level of consciousness, patients with reduced mobility of the arms and/or hands, or patients who are unable to remove the device themselves.

The product shall not be used by patients with a low tidal volume, as the added dead space may cause CO₂ (Carbon dioxide) retention.

1.3 Description of the device

Provox Life™ HMEs are single-use devices for pulmonary rehabilitation. They come in different levels of humidification, breathing resistance and filtration that makes them suitable for different situations.

The different Provox Life™ HMEs are:

Home: when taking it easy, **Go:** when you are out and about, **Energy:** when physically active, **Protect:** when you need protection from bacteria, virus, dust and pollen, **Night:** when sleeping.

Provox Life™ HMEs are compatible with Provox Life™ attachments.

1.4 Technical data

	Pressure Drop at 30 l/min*(Pa)	Pressure Drop at 60 l/min*(Pa)	Pressure Drop at 90 l/min*(Pa)	Moisture loss at VT=1000 ml* (mg/l)
Home	60	200	450	19.5
Go	30	100	230	22.5
Energy	15	50	110	23
Protect**	55	180	320	23
Night	65	210	470	18.5

HME internal volume (Dead space)* max 13 ml

* According to ISO 9360

** Bacterial Filtration Efficiency >98%, Virus Filtration Efficiency >98%. Test method adapted from ASTM F2101.

1.5 WARNINGS

- Do not use other or non-original products as it may cause personal injury or damage the products. Do not alter the product.
- Do not exert unintentional pressure on the lid of the HME. Unintentional or accidental closing of the HME may obstruct breathing.

1.6 PRECAUTIONS

- Always test the function of the HME prior to use. The lid should immediately return to its open position after releasing it.
- Do not disassemble the HME since this will interfere with its proper function.
- The HMEs shall not be used beyond 24 hours, are not reusable and cannot be rinsed in water or any other solution. This will substantially reduce the function of the HME and can increase the risk for infections due to bacterial colonization.
- Do not administer medicated nebulizer treatment over the device since the medication can become deposited in the device.
- Do not use humidifiers or heated humidified oxygen via a mask over the tracheostoma while using the device. The HME will become too wet. If oxygen therapy is required, use only non-heated humidified oxygen.
- Be careful not to exert pressure on Night HME unintentionally. Pressure may reduce the size of the side openings and cause difficulty breathing.

2. Instructions for use

2.1 Operating instruction

For instructions on how to use the products, please see fig. 1-3.

2.2 Device lifetime and disposal

The HME is a single use device and must be replaced at least every 24 hours, or more often if needed. Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

3. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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