
PROVOX[®] LaryTube[™] Life


Rx
ONLY

CE 0413

Figures 1


Figures 2

Ø


↔

8	12.0mm	27/36/55mm
9	13.5mm	27/36/55mm
10	15.0mm	27/36/55mm
12	17.0mm	27/36/55mm


Figures 3


Figures 4a


Figures 4b


Figures 4c


Figures 4d

Figures 4e

Figures 4f

Prescription information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and trademarks

Provox[®] is a registered trademark owned by Atos Medical AB, Sweden. Provox[®] Life[™] and LaryTube[™] are trademarks of Atos Medical AB.

For information about protective rights (e.g. patents), please refer to our webpage www.atosmedical.com/patents.

Manufacturer; Hersteller; Fabrikant; Fabricant; Produttore; Fabricante; Fabricante; Tillverkare; Producent; Produsent; Valmistaja; Framleiðandi; Tootja; Ražotājs; Gamintojas; Výrobce; Gyártó; Výrobca; Proizvajalec; Producent; Producător; Proizvođač; Proizvođač; Κατασκευαστής; Производител; Üretici; მწარმოებელი; Ἀρτοποιητής; İstehsalçı; Производитель; Produzen; Pengilang; 製造業者; 제조사; 製造商; 制造商;

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1. Descriptive information

1.1 Intended use

Provox Life™ LaryTube is a single patient use device intended to provide attachment for Provox Life™ HMEs and accessories after total laryngectomy. For laryngectomized patients with a shrinking tracheostoma it is also used to maintain the tracheostoma for breathing.

1.2 CONTRAINDICATIONS

Provox Life™ LaryTube is not intended to be used by patients that:

- are under any form of mechanical ventilation.
- have damaged tracheal or tracheostoma tissue.

1.3 Description of the device

Provox Life™ LaryTube is a tube made of medical grade silicone rubber. It is delivered single packed and non-sterile. Provox Life™ LaryTube provides attachment for Provox Life™ HMEs and Provox Life™ Shower.

There are 4 models of Provox Life™ LaryTube: (Fig. 1).

Standard model (Fig. 1a) – Made for use with or without a voice prosthesis. Can be attached with a Provox TubeHolder or Provox LaryClips.

Fenestrated model (Fig. 1b) – Made for use with a voice prosthesis. Can be attached with a Provox TubeHolder or Provox LaryClips.

Ring model (Fig. 1c) – Made for use with or without a voice prosthesis. Can only be attached with a Provox Life™ Adhesive.

Fenestrated with Ring model (Fig. 1d) – Made for use with a voice prosthesis. Can only be attached with a Provox Life™ Adhesive.

The Standard and the Ring model can be fenestrated so that air can go through the voice prosthesis. The holes are punched by using Provox Fenestration Punch according to its Instructions for Use. The fenestrations should be in a position that will optimize the airflow through the voice prosthesis.

The different parts of Provox Life™ LaryTube are (Fig. 1e-f):
e) Fenestration holes
f) Ring for adhesive connection

1.4 WARNINGS

- **DO NOT** reuse a used Provox Life™ LaryTube in another patient. Reuse in another patient may cause cross contamination.
- **DO NOT** use lubricating gels if the patient has allergies related to these substances.
- **DO** ensure that the patient has been trained in the use of the device. The patient should demonstrate the ability to understand and consistently follow the Instructions for Use without clinician supervision.

1.5 PRECAUTIONS

- Always use a Provox Life™ LaryTube of an appropriate size. Incorrect size may cause tissue damage, bleeding, irritation, obstruction of the voice prosthesis or may cause the stoma to shrink. HME Cassette may exert pressure on the voice prosthesis if Provox Life™ LaryTube is not sized appropriately.
- **DO NOT** use force during insertion. It may cause tissue damage, bleeding or irritation and accidental dislodgement of the voice prosthesis (Fig.3) and/or XtraFlange if used. If the product is pushed too far into trachea it may obstruct breathing.
- Fenestrate only with Provox Fenestration Punch. This will ensure that the holes are small and smooth. Fenestration holes that are too large may cause dislodgement of the voice prosthesis or XtraFlange. Fenestration holes with sharp edges may cause formation of granulation tissue.
- **DO NOT** use a dirty or contaminated Provox Life™ LaryTube. It may cause infection.
- **DO NOT** clean the device while inside the stoma.
- **DO NOT** use the device until it is completely dry (Fig. 4f). Inhalation of disinfectant fumes can cause severe coughing and irritation of the airways.
- **DO NOT** clean or disinfect by any other method than described as this might cause product damage and patient injury.
- **DO NOT** use a damaged device since this might cause aspiration of small parts or cause formation of granulation tissue.

- Only use water-soluble lubricants. Oil-based lubricants (e.g. Vaseline) may damage the product and should not be used.
- Stop using the device if the stoma gets irritated or starts to bleed.
- Patients with bleeding disorders or patients on anticoagulants should not use the device if it causes bleeding that re-occurs.
- Stop using the device if granulation tissue forms in the tracheostoma.
- Use only original components. Other accessories may cause product damage and/or malfunction.

2. Instructions for use

2.1 Preparation

Choose correct size

- Choose the right diameter and length (Fig. 2).

Prepare for insertion

- Always inspect the product before each use. Do not use if damaged or dirty.
- Make sure that your hands are clean and dry.
- If needed, lightly lubricate the tube with a small amount of the following recommended water-soluble lubricants: KY Jelly[®], SurgiLube[®] and Xylocaine[®] (lidocaine).

2.2 Operating instruction

Insertion (Fig. 3):

Provox Life™ LaryTube Standard and Fenestrated:

Gently insert Provox Life™ LaryTube into the tracheostoma and attach with Provox LaryClip or Provox TubeHolder (Fig. 1).

Provox Life™ LaryTube with Ring:

Make sure that the adhesive is centered properly so that Provox Life™ LaryTube does not press against the sides of the stoma. Gently insert Provox Life™ LaryTube through the adhesive into the tracheostoma (Fig. 1).

NOTE: Insertion and/or removal of Provox Life™ LaryTube may occasionally cause some slight bleeding, irritation or coughing which usually diminishes over time. Contact your clinician if these symptoms persist.

2.3 Cleaning and disinfection

Clean and disinfect Provox Life™ LaryTube before each use. Only use drinking water to clean and rinse the device (Fig. 4a-4f).

1. Clean inside of Provox Life™ LaryTube with Provox TubeBrush or Provox Swab (Fig.4a).
2. The holes of a fenestrated Provox Life™ LaryTube should be cleaned with a Provox Brush (Fig. 4b).
3. Rinse Provox Life™ LaryTube, Provox Brush and Provox TubeBrush in water (Fig. 4c).
4. Place Provox Life™ LaryTube in water-based mild soap for 15 minutes (Fig.4d).

5. Rinse Provox Life™ LaryTube with water (Fig. 4c).
6. Disinfect Provox Life™ LaryTube at least once a day (Fig. 4e) with one of the following methods: Ethanol 70% for 10 minutes or Isopropylalcohol 70% for 10 minutes or Hydrogenperoxide 3% for 60 minutes.
7. Let Provox Life™ LaryTube dry at least for 2 hours or until it looks dry.

2.4 Storage instructions

Store Provox Life™ LaryTube in a clean and dry container at room temperature. Protect from direct sunlight.

2.5 Device lifetime and Disposal

Provox Life™ LaryTube may be used for a maximum period of 6 month. If the device shows any signs of damage it shall be replaced earlier.

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

2.6 Compatible Products

Use only original components. Other products may cause product damage and/or malfunction.

Devices that can be used with Provox Life™ LaryTube are:

- Provox Life™ HMEs (Heat and Moisture Exchangers) for patients breathing through a tracheostoma.
- Provox Life™ Shower is used to prevent water from entering the Provox Life™ LaryTube during showering.
- Provox TubeHolder or Provox LaryClip is used to keep the Provox Life™ LaryTube in place.

- Provox TubeBrush, Provox Swab and Provox Brush
- For Provox Life™ LaryTube with Ring a Provox Life™ Adhesive is used as attachment and keeping the Provox Life™ LaryTube in place.

3. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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© Atos Medical AB, Sweden.