
PROVOX[®] LaryButton[™] Life

Rx
ONLY


CE 0413


42° C
108° F
RT
2° C
36° F


1


2


3


3 a-1


3 a-2


3 b-1


3 b-2


3 b-3


3 b-4


4a


4b


4c


5a


5b


5c


5d


5e

5f

5g

Prescription information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and trademarks

Provox[®] is a registered trademark owned by Atos Medical AB, Sweden. Provax[®] Life[™], LaryButton[™] and LaryClip[™] are trademarks of Atos Medical AB.

For information about protective rights (e.g. patents), please refer to our webpage www.atosmedical.com/patents.

Manufacturer; Hersteller; Fabrikant;
 Fabricant; Produttore; Fabricante; Fabricante;
 Tillverkare; Producent; Produsent; Valmistaja;
 Framleiðandi; Tootja; Ražotājs; Gamintojas;
 Výrobce; Gyártó; Výrobca; Proizvajalec;
 Producent; Producător; Proizvođač; Proizvođač;
 Κατασκευαστής; Производител; Üretici;
 მწარმოებელი; Ἀρτοποιητής; Istehsalçı;
 Производитель; Produsen; Pengilang;
 製造業者; 제조사; 製造商; 制造商;

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1. Descriptive information

1.1 Intended use

Provox Life™ LaryButton is a single patient use self-retaining holder for Provox Life™ HMEs and accessories after total laryngectomy.

For laryngectomized patients with a shrinking tracheostoma it is also used to maintain the tracheostoma for breathing.

1.2 CONTRAINDICATIONS

Provox Life™ LaryButton is not intended to be used by patients that:

- are under any form of mechanical ventilation.
- have damaged tracheostoma tissue.

1.3 Description of the device

Provox Life™ LaryButton is a self-retaining tracheostoma button made of medical grade silicone rubber. The device is delivered single packed, non-sterile and ready to use. The purpose of the device is to create a comfortable and airtight fit between Provox Life™ LaryButton and the tracheostoma, and also to provide attachment for Provox Life™ HMEs and Provox Life™ accessories.

The different parts of Provox Life™ LaryButton are (Fig. 1):

- a) Retention Collar
- b) Shaft
- c) Shield (conical)
- d) HME and Accessory Holder
- e) Wings

1.4 WARNINGS

- Always use Provox Life™ LaryButton of an appropriate size. If Provox Life™ LaryButton is too wide or too long it may cause tissue damage or irritation. If it is too long it may also cause dislodgement of the voice prosthesis or Provox XtraFlange (Fig. 3 a-2). Furthermore, insertion or removal of the HME as well as occluding the stoma by pressing the lid of the HME to speak may exert pressure on the voice prosthesis if Provox Life™ LaryButton is too long.
- Always make sure that the shield (Fig. 1, c) is large enough in relation to the size of the patient's stoma. If the product is pushed too far into the trachea it may obstruct breathing.
- **DO NOT** cut off the wings or otherwise mechanically modify (e.g. fenestrate) Provox Life™ LaryButton since it will change the outer diameter and/or the mechanical stability of the product. This may result in Provox Life™ LaryButton being pushed into the trachea, which may obstruct breathing. A modified Provox Life™ LaryButton may also cause irritation and/or bleeding.
- **DO NOT** reuse Provox Life™ LaryButton between patients. The device is intended for single patient use only. Reuse between patients may cause cross contamination.
- **DO NOT** use force and do not fold the part of Provox Life™ LaryButton that holds the HME during insertion. If the product is pushed too far into trachea it may obstruct breathing and the use of force during insertion may cause bleeding.

- **DO NOT** push the tip of the hemostats or forceps (if used) through Provox Life™ LaryButton (Fig. 3 b-2). If the instrument tip is protruding from the button during insertion, the tissue of the trachea could be harmed and/or the voice prosthesis could become dislodged and subsequently be aspirated or ingested.
- **DO NOT** use lubricating gels if the patient has allergies related to these substances.
- **DO NOT** at any time attempt to clean the device while inside the stoma since this could cause tissue damage and dislodgement of voice prosthesis.
- Make sure that the patient has been trained in the use of the device, and has demonstrated the ability to understand and consistently follow the Instructions for Use without clinician supervision.

1.5 PRECAUTIONS

- **DO NOT** use a dirty or contaminated Provox Life™ LaryButton. Clean and disinfect the device according to the cleaning and disinfection instructions below.
- **DO NOT** use toothed or sharp instruments at any time during insertion since they could damage Provox Life™ LaryButton and cause tissue damage.
- **DO NOT** use a damaged device since this might cause aspiration of small parts or cause formation of granulation tissue.
- **DO NOT** use oil-based lubricants (e.g., Vaseline), they may weaken, damage, or destroy the product. Only use water-soluble lubricants.
- Carefully monitor the tracheostoma tissue during radiation therapy and stop using the device if the stoma gets irritated or starts to bleed.

- Patients with bleeding disorders or patients on anticoagulants should not use the device if it causes re-occurring bleeding.
- Stop using the device if granulation tissue forms in the tracheostoma.

2. Instructions for use

2.1 Preparation

Choose correct size

- Choose the correct diameter and length of Provox Life™ LaryButton.
- To achieve an airtight fit and good retention the diameter of the Retention Collar (Fig. 1, a) of Provox Life™ LaryButton should exceed the diameter of the tracheostoma entrance.
- Ensure that the Retention Collar of Provox Life™ LaryButton fits right behind the edge of the tracheostoma. However if the device is too short it may fall out of the tracheostoma.

CAUTION: During fitting always make sure that Provox Life™ LaryButton does not reach the tracheal flange of the voice prosthesis. Optimal fit and performance of the Provox Life™ LaryButton is achieved if the patient has a so-called “tracheostoma lip” that is situated uniformly around the tracheostoma.

Prepare for insertion

- Always inspect the product before each use. Do not use if damaged or dirty.
- Make sure that your hands are clean and dry.

- If needed, lightly lubricate the tube with a small amount of the following recommended water-soluble lubricants: KY Jelly[®], SurgiLube[®] and Xylocaine[®] (lidocaine).

2.2 Operating instruction

Insertion (Fig. 3):

Alternative 3a:

1. Fold the Retention Collar with your fingers (Fig. 3) and hold it firmly folded.
2. Insert the Retention Collar into the stoma and release it (Fig. 3 a-1 and a-2).

Alternative 3b:

1. Fold the Retention Collar with your fingers (Fig. 3).
2. Use a pair of blunt forceps or a lockable blunt hemostat from the inside to keep it folded (Fig. 3 b-1 and b-2).
CAUTION: DO NOT push the tip of the hemostats or forceps through Provox Life[™] LaryButton (Fig. 3 b-2). If the instrument tip is protruding from the button during insertion, the tissue of the trachea could be harmed and/or the voice prosthesis could become dislodged and subsequently be aspirated or ingested.
3. Insert the Retention Collar into the stoma (Fig. 3 b-3).
4. Release the forceps to unfold the button (Fig. 3 b-4).

Check the correct fit by slightly pulling and/or rotating the button. Gently insert Provox Life[™] HME/accessory into the HME and Accessory Holder.

NOTE: Some patients might experience more coughing after inserting Provox Life[™] LaryButton and when touching the device to speak. This usually diminishes when the user gets used to the device. If discomfort occurs Provox Life[™] LaryButton should be removed.

Device removal

1. Hold Provox Life™ LaryButton in place and gently remove Provox Life™ HME/accessory from the HME and Accessory Holder (Fig. 4a).
2. Hold Provox Life™ LaryButton and gently turn it (about half a turn) to allow the Retention Collar to compress thereby reducing its width (Fig. 4b).
3. Remove the device by gently tilting it out of the stoma (Fig. 4c). Alternatively, the device can be removed by pushing with a finger on the shaft from one direction and tilt the device out.

NOTE: The insertion or removal procedures may occasionally cause some slight bleeding, irritation or coughing. Inform your clinician if the bleeding persists.

2.3 Cleaning and disinfection

Clean and disinfect Provox Life™ LaryButton before each use. Only use drinking water to clean and rinse the device(Fig. 5a-5g).

1. Rinse Provox Life™ LaryButton with water (Fig. 5a).
2. Place Provox Life™ LaryButton in water-based mild soap for 15 minutes (Fig. 5b).
3. Clean inside of Provox Life™ LaryButton with Provox TubeBrush or Provox Swab and rinse with water (Fig. 5c + 5d).
4. Disinfect Provox Life™ LaryButton at least once a day with one of the following methods: Ethanol 70% for 10 minutes or Isopropylalcohol 70% for 10 minutes or Hydrogenperoxide 3% for 60 minutes (Fig. 5e).
5. Let Provox Life™ LaryButton dry for at least 2 hours or until it looks dry (Fig. 5f).

6. Rinse Provox TubeBrush or Provox Swab in water after use (Fig. 5g).

2.4 Storage instructions

Store Provox Life™ LaryButton in a clean and dry container at room temperature. Protect from direct sunlight.

2.5 Device lifetime and Disposal

Provox Life™ LaryButton may be used for a maximum period of 6 month. If the device shows any signs of damage it shall be replaced earlier.

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

2.6 Compatible Products

Use only original components. Other products may cause product damage and/or malfunction.

Devices that can be used with Provox Life™ LaryButton are:

- Provox Life™ HMEs (Heat and Moisture Exchangers) for patients breathing through a tracheostoma.
- Provox Life™ Shower is used to prevent water from entering Provox Life™ LaryButton during showering.
- Provox Life™ BasePlate Adaptor, to attach products with an ISO 15 mm standard connector to Provox Life™ LaryButton.
- Provox TubeHolder or Provox LaryClip can be used to keep Provox Life™ LaryButton in place.

3. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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© Atos Medical AB, Sweden.